

SkyAngel *Classic*™

FRIENDS User Guide



Charging:

Use the **Magnetic Charging Cable** to charge the unit for **3 hours** initially to allow the unit battery to reach 100%.

If the unit is turned ON, when you attach the **Magnetic Charging Cable** the unit will speak "**Charger is connected**".

If the unit is turned ON, when the **Magnetic Charging Cable** is removed the unit will speak "**Charger is removed**".

The **front light** will turn **Red** when charging. The front light will turn **Green** when fully charged.

Press the **Red SOS button** once quickly and you will hear the unit speak the **Battery Charger Percentage**. The unit will speak "**Battery is xx**" xx represents the percentage of charge.

The **battery** can last up to **4 days** at full charge.

For **the best results**, charge about **30 minutes every day** to ensure the unit is always ready to use.

Power ON – Power OFF:

To turn the unit ON.

If the unit is off, **press and hold** the **Side Button** or the **Red SOS Button** for **15 seconds** to power the unit on. You will hear the unit play a chime sound.

As the unit is turning on you will see the **Front Light** turn **Red**, then **Green**, then **Blue**.

When the unit has completely turned on and fully connected to a cell tower there will be no lights on the front of the unit. When you quickly press the SOS you will hear the time and battery percentage. This confirms the unit is connected and ready.

To turn the unit Off.

From the App Main Screen select **Remote Restart** or **Shutdown The device**.

SOS Alert Functions:

This unit has THREE methods to activate the SOS function and call 911:

Method #1 Red SOS Button:

When the unit is powered on, press and hold the **Red SOS Button** or the **Side Button** for 3 seconds. You will feel the unit vibrate and you will hear the unit speak “**Sos starting**”. The unit will then call your Three emergency contacts in the order they were entered in the APP. The unit will call the first number and wait approx. 15 seconds for an answer, if there is no answer it will hangup and call the second number, and so on. It will repeat this calling cycle two times. Once your contact answers you can talk to them through your unit and the unit will not call the next number in sequence.

Method #2 Selective Calling Red SOS Button:

Three contacts can be called in order they are entered in the APP.

If you quickly **press the Red SOS Button or the Side Button Twice** it will skip the first number and directly call the Second number in your contact sequence.

If you quickly **press the Red SOS Button or the Side Button Three times** it will skip the first two numbers and directly call the Third number in your contact sequence.

**Hint: This allows you to call the second and third numbers directly without waiting for the unit to call the numbers before that.*

**How to use idea: This would be helpful if the third number is 911 and you would want to skip the contacts and call 911 first.*

Method #3 Automatic Fall Detection:

If the unit **detects an injury causing fall of 2 feet or more above the ground**, you will hear the **unit speak “Sos starting”** then the unit will automatically call your contacts. You do not need to press any buttons. The call will be placed automatically. Once your contact answers you can talk to them through the unit.

To Cancel an Sos Call:

Press and Release the Side Button.

****Please allow the unit up to 2 minutes to connect a call if the cellular network is busy****

Answer incoming Calls:

This unit has the ability answer incoming calls. When the unit rings press the **Green** button to answer. Others would simply call the units phone number the same way they would call any number. The telephone number of your unit will be provided during the activation process.

Send/Receive Short Voice Messages:

Sending a voice message from the unit: To send a short voice message (up to 15 seconds) to your APP simply press and HOLD the **Green** button for 5 seconds or until you hear “Voice recording in progress”, keep HOLDING the **Green** button while speaking your message. When done release the **Green** button. You will hear a confirmation sound.

Listening to received voice messages on the unit: To listen to messages that have been sent to the unit press the **Green** button 1 time quickly. Messages are deleted manually within the APP. Voice messages cannot be deleted on the unit.

Accessing voice messages from the APP: To listen to a message, access the APP then click CHAT. You will see the messages there. Click the **Green** message bubble to hear it. Press and hold the **Green** message bubble to Delete it.

To record a voice message, access the APP and click CHAT. At the bottom of the screen you will see ‘START RECORDING’, press and hold the button while recording your

message. Let go when done and the unit will play a tune to notify the user there is a message waiting for them on the unit.

**Hint: Giving friends and family access to the app will allow them to send and receive voice messages from the unit, track the location of the unit and see health data. You can provide View only access or Full access. Full access allows for changing of settings. View only access allows for viewing only.*

Green Button Health Features:

Keep in mind to achieve more accurate results the American Medical Association recommends:

Sit in a supported chair with back straight, feet flat on the floor, and legs uncrossed for at least 5 minutes before the measurement. Your arm should be supported on a flat surface at heart level. You may hold your watch hand over your chest with your other hand.

Pulse:

Press the **Green** Button **Twice** Quickly:

You will hear **"Testing pulse please wait a moment"**, in about 1 minute and 30 seconds you will hear the result as Example: **"Pulse 60 High Pressure 120 Low Pressure 73"**

Oxygen or SP02:

Press the **Green** Button **3 Times** Quickly:

You will hear **"Testing Sp02 please wait a moment"**, in about 1 minute you will hear the result as Example: **"Sp02 97%"**

Temperature:

Press the **Green** Button **4 Times** Quickly:

You will hear **"Testing body temperature please wait a moment"**, in about 30 seconds you will hear the result as Example: **"The body temperature is 87 point 3"**

Cancel a Health Measurement during the testing:

Press the **Side Button** quickly and release to cancel the measurement. The unit will speak **"Cancel Measurement"**.

Other:

Programming the Contacts:

3 emergency contacts can be programmed by accessing the Smartphone APP associated with the unit.

Smartphone APP features:

Track The Unit on a Map – See the current location of the unit on a map.

Create a Phonebook – can be used to allow certain numbers to call you

Health Monitoring – all of the health features such as Steps, Sleep, Heart Rate, Blood Pressure, Blood Oxygen, Body Temperature can be reviewed over time to see trends.

Automatic Monitoring – your vital signs can be automatically monitored at regular intervals such as each hour so you can see a regular history of the results.

The following setting can be adjusted using the APP; Fall Detection Sensitivity, Time Zone, Language, Family Members, Text Messages, Call Log, Block Unknown Callers, Set Time Format, Adjust Volume, Schedule Automatic Reboot, 3 Alarm Reminders, Auto Answer, Find Lost Unit by making a sound on the unit, Set Vibration and Ring, GPS update frequency, Lock if Lost, Set Do Not Disturb Times.

Battery Level:

Press the **Red Sos Button** once quickly and the unit will speak the current time and then the battery level. Example: **“The time is 7:15pm Battery 90”**. 90 represents the percentage of battery level.

Waterproof:

You can shower or take a bath with this unit for up to 30 minutes. After 30 minutes of submersion please wait 24 hours before re-submersing. It is IP67 certified waterproof.

Wrist Band is replaceable:

The wrist band size is 19 millimeters or $\frac{3}{4}$ inch. You can easily replace the band with another strap design by removing the pins and replacing the band.

Cellular:

This unit uses all available **3G** and **4G** cellular towers and carriers.

Size and Weight:

Length of unit from top to bottom: (Including band) 9.5 inches

Width unit face: **1 ¼** inches - **Top to bottom** of face: **1 7/8** inches

Thickness of unit: **5/8** inch - **Weight: 1.6** ounces

Battery replacement:

The battery is 600mah. The battery has a 1-year warranty from the factory. If the battery needs replacing outside the warranty, contact a certified battery replacement specialist or contact ATSCares.com.

Device Tag:

A device tag as seen below is attached to the watch. Please cut the clear plastic loop to remove. Please retain the tag in the event you would want to upgrade the watch in the future you will need the information on the tag.

Example Tag



Warranty:

This product is protected by a one-year manufacturer defect warranty. Please contact AtsCares.com to activate warranty coverage. Items covered under warranty are to be shipped to the warranty center paid by customer. The Item will be repaired or replaced under warranty. Warranty center will pay for return shipping in the U.S. Coverage outside the U.S., all shipping costs will be paid by customer.



SkyAngelClassic.com

Scan this QR code with your camera for information on your SkyAngelClassic