

AMPLIFIED EMERGENCY MEDICAL CONNECT

**With emergency remote calling - - - and
NO MONTHLY FEES or CONTRACTS - - - EVER!**

MEDICAL ALERT SYSTEM with 2 waterproof wrist or necklace pendants



Catalog #: ats600gm

This combination **Amplified BIG button Phone** with **Caller ID, Speakerphone** and **Emergency Connect** provides peace of mind for those who live alone and their families.

Includes: 2 waterproof pendants that can be worn on the wrist or around the neck and has an amazing list of features with 1 Year Warranty

ONLY \$159.99 + sales tax / shipping & handling

*Prices subject to change without notice
ships in 2 business days

The bonus of an amplified phone WITH addition of the emergency connect feature. The Amplified Emergency Medical Connect Phone acts as an automatic dialer in emergency situations. It comes with 2 small transmitter remotes that can be worn on wrist like a watch, or attached to a lanyard around the neck (included).

In the event of an emergency, the wearer presses the red emergency button on the transmitter OR on the phone. Once the emergency button is pressed, the phone numbers of up to six programmed emergency contacts are dialed. These are normally family members, friends or neighbors who would be able to respond to the emergency call. If there is no response after 30 seconds or the line is busy, the phone automatically dials the next programmed number. It will cycle through numbers twice. When the phone reaches a live person, it will play your pre-recorded emergency message. The other person will dial a number (0-9) to confirm that the emergency message has been received and to deactivate the system from dialing the next emergency contact. The system then turns on its speakerphone so the other person is able to speak and listen to the other user, known as Remote Audio Monitoring.

This system also has a remote call feature that allows a caller to call the phone and after 10 rings enter a user password that will **remotely answer the phone in 2 way speakerphone mode**. This allows you to call and check on someone and speak to them through the speakerphone without them having to touch the phone. This is a great way to call and make sure they are ok without the person you were calling having to do anything to answer the phone, totally hands free. Easy to install and use - - phone uses standard phone line - - transmitter pendants are battery-powered, allows you to remain independent - - provides added security without monthly monitoring fees from security companies

FEATURES: 50+ dB amplification - - Remotes that operate up to 100 feet from the base unit - - Up to 6 remote-dial contact numbers - Emergency dial button - - Number displays on screen during dialing - Strobe light ring flasher - - Super-loud 95dB ringer (adjustable volume) - - Speakerphone - - Auto on-hook after 8 minutes - - Call cancel on pendants and base - - Caller ID (99 name and number storage) - - Phonebook (99 name and number storage) - - Go to next number every 30 seconds - - Remote audio monitoring (ability for caller to force the phone to answer in speakerphone mode) - - Adjustable amplification with volume control slide bar, amplification boost button and situational boost button - - Adjustable tone control with full amplification in high and low tones - - Amplified 3.5mm jack /

Amplified 2.5mm hands-free jack - - Hearing aid t-coil compatible - - 4 Ringer tones - - 6 Memory buttons and 3 priority memory buttons - - Up to 12dB speech output amplification - - Bed shaker (optional - sold separately) - - Desk or wall mount - - Adjustable angle display - - Tone or Pulse dialing - - Remote requires one 12V battery (replaceable) - - Requires a standard 120V electrical outlet - - *One Year Limited Warranty*

INCLUDED IN PACKAGE:

2 Wireless Waterproof Pendant Transmitters

3 Wrist Straps - - 2 Necklace Lanyards

2 Telephone Cord 1-Long - - 1 Short

1 Wall mount bracket

1 Wall Power Adapter



ORDER FORM

Assistive Technology Services

"Helping You Live Independently"

8023 Franklin Road
 Murfreesboro, TN, 37128
 (615) 562-0043

Billing Address Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing []	Ship To Address Name _____ Address _____ City _____ State _____ Phone Number _____
---	---	---

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	Ats600gm	2 Wrist/Necklace Pendant Medical Alert System US Shipping	\$183.77		
	Ats600gm	2 Wrist/Necklace Pendant Medical Alert System Canada Shipping	\$200.78		
	Ats600gm	2 Wrist/Necklace Pendant Medical Alert System Shipping to Other	\$214.87		
				Subtotal	
		Only Applicable If Shipping To Tennessee		TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice		Shipping	
				Total	

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [] MC [] AmEx [] Disc []				
Name On Card	Email Address			

If Paying By Check Please Make Checks Payable To: Assistive Technology Services

"May You Know The Glory Of God And May He Bless You"