

# Fall Guard Alert System



*The Fall Guard Alert System will alert you if your loved one leaves their bed or chair.*

*It has a tamper-proof monitor that can only be deactivated with a key. It also has a voice prompt.*

*Great for small facilities, caregivers or at home. Anywhere you need the help of the Fall Guard Alert System.*

**ONLY \$149.62 With Chair Pad** US + Shipping/Handling

**ONLY \$164.09 With Bed Pad** US + Shipping/Handling

\*Prices are Subject to Change Without Notice

**ITEM # ats4520**

## **Do you worry about your loved one falling or getting out of bed while you sleep?**

The Fall Guard Alert System has a disposable Pad that triggers the alarm if weight is removed from it.

The Fall Guard Alert Systems only has two parts.

### **The Alert Monitor with key and Pressure Pad (Disposable)**

The Alert Monitor can be powered by a 9 volt battery (not included) or it can be plugged in the wall with an AC adapter (included).

If you choose the Caregiver's Key option, the Alert Monitor can only be deactivated with the Key. This way your loved one cannot turn it off.

The Voice Prompt can be set to alert your loved one to "Please lie back down" or "Please sit back down" before the alarm is triggered.

This Voice Prompt is usually enough to get your loved one to lay or sit back down.

It has an Automatic Reset that prevents you from forgetting to turn it back on after the alarm has been silenced.

The controls for the Alert Monitor are made to be difficult to get to. This way your loved one cannot get to them and change the settings or turn it off easily.

This Fall Guard Alert System is a hospital grade product. They are as good at your home as they are in a nursing home or hospital. This is a product you will love and trust. The professionals do!

The pads are disposable and should only be used for 90 days. (Extra bed pads and chair pads are available.)

### **Details**

Dimensions Bed Pad: 5" x 30"

Chair Pad: 7" x 15"

Alert Monitor: 6" L x 3½" W x 1¼" H

Product Weight 0.95 lb



**Bed Pad**

**Chair Pad**





# ORDER FORM

**Assistive Technology Services**  
*"Helping You Live Independently"*

8023 Franklin Road  
 Murfreesboro, TN, 37128  
 (615) 562-0043  
 Fax (831) 308-9671

<b>Billing Address</b> Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing [ ]	<b>Ship To Address</b> Name _____ Address _____ City _____ State _____ Phone Number _____
---	---	---

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	Shipped to US	Shipped to Canada	Shipped to Other	UNIT TOTAL
	Ats4520	Fall Guard Alert System with Chair Pad	\$156.51	\$176.51	\$196.51	
	Ats4520	Fall Guard Alert System with Bed Pad	\$170.98	\$190.98	\$210.98	
					<b>Subtotal</b>	
		Only Applicable If Shipping To Tennessee			TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice			Shipping	
					<b>Total</b>	

**If Paying By Check Please Make Checks Payable To: Assistive Technology Services**

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [ ] MC [ ] AmEx [ ] Disc [ ]				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"