

# Home Support Pole-Grab Bar



**Need support getting out of bed, in the shower, in the bathroom, down steps?**

**This will give you something to hold on to while getting up, down and around.**

**This floor to ceiling support pole is for you. Compare our prices to those that are \$300 or more.**

**ONLY \$99.00** US + Shipping/Handling

\*Prices are Subject to Change Without Notice

**ITEM # ats4391tp**

Galvanized steel tubing is designed to support you easily and with portability.

NOT permanently attached. Easy to move around.

Can be used anywhere in the home.

The ductile iron rack, handle and cast aluminum housing are proven performers for durability and support.

The durable molded neoprene swivel feet provide excellent grip in wet environments like the shower or outdoors.

Pivoting feet to provide more gripping power.

Adjustable Range of 89" - 104" (7 feet 4 inches to 8 feet 6 inches)



**Assistive Technology Services - [www.ATS-TN.com](http://www.ATS-TN.com) - (615) 562-0043**



# ORDER FORM

**Assistive Technology Services**  
*"Helping You Live Independently"*

8023 Franklin Road  
 Murfreesboro, TN, 37128  
 (615) 562-0043  
 Fax (831) 308-9671

<b>Billing Address</b>  Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing [ ]	<b>Ship To Address</b>  Name _____ Address _____ City _____ State _____ Phone Number _____
---	---	---

CUSTOMER EMAIL ADDRESS		PAYMENT TERMS		DATE
		Due on receipt		

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	Ats4391tp	Home Support Pole-Grab Bar Shipped to US	\$99.00	\$17.99	
	Ats4391tp	Home Support Pole-Grab Bar Shipped to Canada	\$99.00	\$37.99	
	Ats4391tp	Home Support Pole-Grab Bar Shipped to Other	\$99.00	\$57.99	
				Subtotal	
		Only Applicable If Shipping To Tennessee		TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice		Shipping	
				<b>Total</b>	

If Paying By Check Please Make Checks Payable To: Assistive Technology Services

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [ ] MC [ ] AmEx [ ] Disc [ ]				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"