

# Lighted Cane



*The Lighted Cane has a bright white light that illuminates your surroundings. Push the button on the handle and the Xenon bulb inside will shed a bright light through the clear bottom half of the cane.*

*Offset handle.*

*Perfect for negotiating steps or illuminating the path in a dark garage.*

**ONLY \$43.99** US + Shipping/Handling

\*Prices are Subject to Change Without Notice

**ITEM # ats4531**

Do you need some extra light for a midnight trip to the bathroom? How about walking to the door at night? If you use this cane you won't need a nightlight or a flashlight! This cane is elegant by day and provides the light you need to see where you are going when it is dark!

The top half of the lighted cane has a shiny black finish and an offset handle with a foam cover for a comfortable, stable grip.

Push the button on the end of the handle and the Xenon bulb inside will shed a bright light through the clear bottom half of the cane. The light covers a 3' to 4.5' diameter, perfect for negotiating steps or illuminating the path in a dark garage.

## **Specifications:**

Dimensions Stem: 36" L x 3/4 " diameter

Handle: 1 1/2" diameter

Weight Capacity 250 lb.

Product Weight 1.5 lb.

Warranty One Year Manufacturer's Warranty

Required 3 "AA" Batteries (not included)





# ORDER FORM

**Assistive Technology Services**  
*"Helping You Live Independently"*

8023 Franklin Road  
 Murfreesboro, TN, 37128  
 (615) 562-0043  
 Fax (831) 308-9671

<b>Billing Address</b> Name _____ Address _____ City _____ State _____ Phone Number _____	Check if same as Billing [ ]	<b>Ship To Address</b> Name _____ Address _____ City _____ State _____ Phone Number _____
---	---	---

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	Ats4531	Lighted Cane Shipped to US	\$43.99	\$8.98	
	Ats4531	Lighted Cane Shipped to Canada	\$43.99	\$28.98	
	Ats4531	Lighted Cane Shipped to Other	\$43.99	\$48.98	
				<b>Subtotal</b>	
		Only Applicable If Shipping To Tennessee		TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice		<b>Shipping</b>	
				<b>Total</b>	

**If Paying By Check Please Make Checks Payable To: Assistive Technology Services**

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [ ] MC [ ] AmEx [ ] Disc [ ]				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"