

Take your “Emergency Alert System” with you wherever you go! 24 /7 anyWHERE – anyTIME protection with

Stay protected
24 / 7

SkyAngelGPS™



You don't have to pay each month for a medical alert system you can only use at home? You can now have the same safety, comfort and convenience PLUS more, for the same price ... with a GPS medical alert from SkyAngelGPS™.

Now you can take your protection anywhere you go – take it to the lake, to town, have in your car, on vacation, shopping – anywhere you go! This is the great advantage of a GPS medical alert system from SkyAngelGPS™.

Free yourself from your home and keep the comfort knowing you have an *Angel* watching over you! SkyAngelGPS™ is always with you.

“What makes SkyAngelGPS™ different?”

SkyAngelGPS™ is similar to an “OnStar for Senior Citizens”. This premier GPS medical alert system provides freedom with security. With the touch of a button, you can talk with help, anytime, from virtually anywhere. Other products like a traditional Personal Emergency Response System (PERS) only allows you to call for help within 150 to 600 feet away of a base station. Unlike a PERS, SkyAngelGPS™ works both indoors and outdoors to ensure anyone has peace of mind knowing that help is always available. In addition, the SkyAngelGPS™ has GPS enabling the emergency call center to locate you, even if you're lost or can't describe your location.

“Is it complicated to operate?”

NO! This unit has a one-button instant contact feature with a 2-way hands-free speakerphone also! Press 1 button – that's it – and your SkyAngelGPS™ is activated and the 2-way speakerphone allows instant communication with the call center help specialist. Simple as that!

“Why not just use a cell phone instead of the SkyAngelGPS™?”

The expectation is in an emergency, you will have the cell phone fully charged and with you, be able to open it up and dial 911 AND be able to communicate your name, exact location and emergency. The advantage is you may already have a cell phone, however, it is a poor substitute for a medical alert system and provides a false sense of security.

There are 4 basic problems calling 911 from a cell phone:

- 1- They don't know WHO you are: The caller ID from your cell phone does not identify who you are, only your phone number and possibly a general metropolitan area where your cell service originates from. They will not know your name, home address, medical condition or who to call in the event of an emergency.
- 2- They don't know WHERE you are: 911 responders don't immediately know where you are – the closest they know is the location of the cell service tower you used when placing your call. If you lose connection or drop your phone before you have given your location, they may not be able to find you.
- 3- All cell phones DO NOT send GPS coordinates. Your phone may have a GPS feature on it, but that information may not be sent to 911.
- 4- When an accident, emergency or medical situation occurs, finding the buttons and dialing an emergency number could be difficult under stressful situations. Having a 1-button operation solves this problem.

“Can I try it and return it if I am not satisfied?”

Yes. There is a 30-day return policy. You can try it for 30 days and return it for a refund if you are not satisfied. All shipping costs are the responsibility of the customer including our original shipping/handling costs incurred to ship the SkyAngel/GPS™ to you. These shipping costs will not be refunded with a return.

“Can I use the SkyAngel/GPS™ if I have a pacemaker?”

Yes. The FDA has determined they “would not seem to pose a significant health problem ... to pacemaker wearers.”

Read the report at:

<http://www.fda.gov/Radiation-EmittingProducts/RadiationEmittingProductsandProcedures/HomeBusinessandEntertainment/CellPhones/ucm116311.htm>

“How does it work?”

SkyAngel/GPS™ is simple to use. Press the emergency button and it automatically places a call to our around-the-clock professional responders. You don’t need to remember numbers or have to dial a number --- it happens at the press of a button. Emergency personnel stay online with you until appropriate help arrives at your location. Your family, caretaker and doctor are immediately notified.

“How do I set it up?”

Easy! The SkyAngel/GPS™ device does not require any setup on your part. We perform quality checks and take great care in preparing it for use. The only thing you need to do with the device is make sure it’s charged and turned on.

“How do I know if SkyAngel/GPS™ will work in my area?”

SkyAngel/GPS™ uses T-Mobile cellular network – visit <http://www.t-mobile.com/coverage/pcc.aspx> and enlarge the coverage map to find your area.

FEATURES:

- Works anywhere there is cell phone service
- 1-Button call for assistance
- Instant 2-way speakerphone communication
- GPS tracking allowing emergency help to know exactly where you are
- Works 5-12 days, between charges, upon a full charge
- NO Contracts – NO Commitments
- FREEDOM to leave your home and carry the *Angel* with you!
- One Year FULL WARRANTY on all equipment

SPECIFICATIONS	
Weight	3.1 ounces
Operation Time	5-12 days depending on continuous tracking being enabled and check-in time
Internal Battery	Rechargeable and removable - 850mAh battery (3.7V)
Charging	DC 4.2-5.5V/400mA (mini USB port)
GSM Module	GSM 850/1800/1900Mhz
GPS Chipset	Assisted GPS with LocalTouch™ Technology for Enhanced Accuracy

“How can I get one?”

- * Complete an Order Form and mail with check or credit card information;
- * Visit our website at <http://www.skyangelgps.com/> and pay securely online;
- * Call our support line at 615-562-0043 Mon-Fri 8am – 5pm CST



ORDER FORM



www.SkyAngelGPS.com

Assistive Technology Services

"Helping You Live Independently"

8023 Franklin Road
Murfreesboro TN 37128
Phone: (615) 562-0043

Billing Address Name _____ Address _____ City _____ State _____ Zip _____ Phone _____	Check if Ship to is same as Billing <input type="checkbox"/>	Ship To Address Name _____ Address _____ City _____ State _____ Zip _____ Phone _____
--	---	--

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	atsSkyAngelGPS	Monthly Plan	99.00	included	
		Includes setup charge, shipping, activation, SkyAngelGPS™ Alert, rechargeable battery, charging cord, power adaptor, lanyard & first month	+ 39.99 /mo		
	atsSkyAngelGPS	Annual Plan (save over \$100)	449.99 /yr	included	
		Includes setup charge, shipping, activation, SkyAngelGPS™ Alert, rechargeable battery, charging cord, power adaptor, lanyard & one year			
	NOTE:	After purchase, you will receive a contract which must be executed prior to activation			
				Subtotal	
		Only applicable if shipping to Tennessee		TN Sales Tax 9.75%	
		*Prices are subject to change without notice		Shipping	
				Total	

If paying by check - make checks payable to: Assistive Technology Services

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code (back)
Visa <input type="checkbox"/> MC <input type="checkbox"/> AmEx <input type="checkbox"/> Disc <input type="checkbox"/>				
Name On Card:		Billing Email Address:		

"May You Know The Glory Of God And May He Bless You"