


Super Pole With Bar



Add stability with a horizontal bar to hold on to.

The Super Pole with bar provides the additional supporting feature of a horizontal bar that moves.

See video for more details

ONLY \$259.99 US + Shipping/Handling

*Prices are Subject to Change Without notice

ITEM # ats22015

The superpole with bar gives you the **extra support** to grab something when getting in our out of a tub.

The **horizontal bar** will swivel 360' around the pole and can 'lock' in place every 45 degrees for stability.

The horizontal bar can also be adjusted from 28 inches to 58 inches above the ground. This allows to set the height of the bar to exactly the right height that is needed.

The superpole is a great solution when you do not have a surface to install a grab bar in the bathroom.

The superpole is height adjustable from 93 to 99 inches. The pole is 1.5 inches in diameter.

The pole is **held solid in place** by an easy to use 'jack-screw' between the floor and the ceiling.


Rubber pads keep the pole from scarring up the floor or ceiling.

Easy to move if you need to. No drilling, no screws, not permanently mounted. **Take it with you !**

Can hold up to **300 pounds**.

Returns are subject to a 20% restocking charge due to size and re-packaging of item.

Comes with a 1 year limited warranty.



Assistive Technology Services - www.ATS-TN.com - (615) 562-0043



ORDER FORM

Assistive Technology Services
"Helping You Live Independently"

8023 Franklin Road
 Murfreesboro, TN, 37128
 (615) 562-0043
 Fax (831) 308-9671

Billing Address Name _____ Address _____ City _____ State _____ Phone Number _____		Check if same as Billing [<input type="checkbox"/>]	Ship To Address Name _____ Address _____ City _____ State _____ Phone Number _____
---	--	--	---

CUSTOMER EMAIL ADDRESS	PAYMENT TERMS	DATE
	Due on receipt	

QTY	ITEM #	DESCRIPTION	UNIT PRICE	SHIPPING	UNIT TOTAL
	Ats22015	Super Pole with Bar Shipped to US	\$259.99	\$18.67	
	Ats22015	Super Pole with Bar Shipped to Canada	\$259.99	\$38.67	
	Ats22015	Super Pole with Bar Shipped to Other	\$259.99	\$58.67	
				Subtotal	
		Only Applicable If Shipping To Tennessee		TN Sales Tax 9.75%	
		*Prices are Subject to Change without Notice		Shipping	
				Total	

If Paying By Check Please Make Checks Payable To: Assistive Technology Services

Credit Card Type	Credit Card Number	Expiration Month	Expiration Year	Security Code From Back Of Card
Visa [<input type="checkbox"/>] MC [<input type="checkbox"/>] AmEx [<input type="checkbox"/>] Disc [<input type="checkbox"/>]				
Name On Card	Email Address			

"May You Know The Glory Of God And May He Bless You"