

Help is always within reach.

A two-way voice medical alert pendant with automatic fall detection - designed for simple, dependable help at home.

NO MONTHLY FEES

No monitoring subscription required.



Simple safety. Clear communication.

2W

Two-Way Voice on the Pendant

Talk and listen directly through the wearable pendant - no need to reach a separate phone.

FD

Automatic Fall Detection

If a fall is detected, the pendant can begin the help-call sequence automatically, even if you cannot press the button.

SOS

SOS: 3 People + 911

The SOS sequence can call up to three selected contacts, then 911 as the emergency number - up to four numbers total.

H2O

Waterproof Pendant

IPX65 rated for shower use. Keep the USB cover secure and do not submerge the pendant.

900

Long-Range Protection

Designed to operate up to 900 feet from the base station, depending on the home environment.

1

One-Button Help

A large tactile button makes it easy to start an emergency call without navigating menus or small controls.

TEL

Home phone line required for standard setup.

Works with compatible traditional, cable, or Internet phone service. No landline? Ask us about alternative solutions.

How Touch n Talk FD gets help

Whether you press SOS or the pendant detects a fall, the system follows a simple call sequence.

Emergency call sequence

1

Wear the pendant

Keep the two-way voice pendant with you around the home and wear it in the shower.



2

Press SOS or fall is detected

Hold the call button to request help, or let automatic fall detection trigger the sequence.



3

Calls up to 3 trusted people

Program family, friends, neighbors, caregivers, or other chosen help contacts.



4

Then calls 911

If help contacts do not connect, the emergency number can be set to 911 - up to 4 numbers total.

Two-way voice lets the person receiving the call speak directly with the wearer through the pendant.

Everything needed to get started

The standard system includes the core components for wearing, charging, and connecting the pendant to your home phone service.



Also included: USB charging cable. The pendant is rechargeable and gives spoken voice prompts to guide the user.

Need more coverage? The system can support multiple compatible voice pendants. Ask ATS about additional pendants and optional accessories.

Designed for everyday independence

No monthly fees

Uses your home phone service instead of a monitoring subscription.

Two-way speakerphone

Speak and listen through the pendant during a call.

Voice prompts

Spoken guidance helps explain what the device is doing.

Up to 100 contacts

Beyond emergency contacts, the pendant can also function as a simple cordless phone.

Shower use

IPX65 rated. Do not submerge; secure the USB port cover.

Alternative solutions

No landline available? Contact ATS to discuss connection options.

Questions? 615-562-0043
atscares.com